**Schüttelwasserbäder
Shaking Water Baths**

1083 / 1086

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**Operating
Instructions**

**1083 / 1086 Schüttelwasserbad
Shaking Water Bath**

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
1. Before installation

The information given in the present manual must **by all means** be carefully read and observed. Only then can a perfect functioning of the Shaking Water Bath be guaranteed.

**Vital information within the manual are emphasized in bold letters.
Safety precautions are additionally marked with the following symbols**

Warning of hot surfaces

**Warning of dangerous
electrical voltage**

A free of charge guarantee repair cannot be granted for defects due to improper installation or handling.

2. Location of the Shaking Water Bath

Place on solid, even and level surfaces only. **Do not use outside buildings.** The unit is **not** suitable for use in explosion endangered surroundings, eg during anaesthesia with inflammable gas or steam types!

3. Voltage

Main switch and starter switch must be in **OFF** position.

Mains voltage and voltage stated on the name plate at the back of the unit must be identical.

If this applies, unit can be connected.

The Shaking Water Bath must only be connected to a correctly installed shock proof power socket.

4. Installation

The Shaking Water Bath must only be operated with tap water. Other media, such as eg oil, acid or agents against bacteria, etc (chlorine or copper sulphate) can lead to damages of the bath interior, the screwing or the heating element.

Before operation, the bath interior must be filled with water. Do make sure that the water outlet tap at the back of the unit is closed (handle must have a 90° angle from tap position). Open the lid and fill the bath. The water level should be kept between the "**min**" and "**max**" markings.


5. Special Accessory Level Regulator

With the adjustable level regulator at the back of the unit, occurring water loss (through evaporation) is compensated and thus water level is kept constant. The level regulator also allows operation of the bath below room ambient, minimum of approx. 3°C above tap water temperature. Before operation, the level regulator must be connected to water tap.

The upper hose connection of the switch is the water supply. Connect this to the tap by means of a laboratory hose (inside diameter max. 9 mm). The lower hose connection is for water outlet and must be connected to a lower placed drain (outlet descendent) by a hose as described above. Please ensure a free flow for this outlet without the possibility of backflow!

All hose connections must be secured with hose clamps.

The desired water level can be adjusted by the outlet tube of the level regulator. For this purpose please loosen the outlet tube screwing with a spanner GW 27. The water level can now be set by pulling out or pushing in the tube. Tighten screwing again. Now carefully open water tap and fill Shaking Water Bath with water as described in article 4.

6. Starting Operation

Switch on main switch. The green pilot lamp will glow.

7. Temperature Regulation

After switching on the Shaking Water Bath, the display shows the actual inside temperature. By pressing key "°C", the display will show the last set and registered temperature that will now be used again automatically.

The temperature range is between approx. 5°C above ambient to 99.9°C, temperature constancy (temporal) approx. +/- 0.1°C. The set temperature can be changed with keys "°C", "+" and/or "-". By pressing key "°C", the display switches from actual to set temperature. By pressing key "°C" and either "+" or "-" simultaneously (two finger operation), the desired temperature can be set in 0.1°C steps. The counting speed increases after a short time. The newly set temperature is saved as soon as key "°C" is released. The display will now show the actual temperature again.

The yellow pilot lamp next to the temperature display shows that heating is in operation.

CAUTION:
**Hot surfaces when temperatures rise
above 60°C
Do not touch!
Danger of burns!**

8. Cooling Coil (only Shaking Water Bath 1086)

The standard cooling coil, suitable for connection to water taps or usual external cooling appliances, extends the unit's temperature range (lowest possible operation temperature: +10°C). The cooling coil connection is situated at the back of the unit. Hoses of an internal diameter of approx. 6 mm can be pushed onto the cooling coil connection tube – **hoses must be secured with hose clamps.**

9. Error Display

9.1. The temperature regulator disposes of electronic monitoring with error screening and concurrent over-temperature cut-out.

The **Temperature Sensor** of the regulator is constantly being screened for short circuits. In case of fault, display will show: "**E1L**" for a short circuit in the Temperature Sensor or "**E1H**" for interruptions. The microprocessor controlled temperature regulator will switch off the heating circuit.

After elimination of error, unit can be restarted.

The unit also disposes of an **electronic over-temperature cut-out** for the protection of test substances, dependent on the set temperature. This device will switch off heating when the set temperature is exceeded by 4°C. The electronic monitoring will induce the information "**E2H**" to be shown on the display.

To restart heating, the unit must be switched off and then on again at the main switch. Should the information "E2H" remain, error must be attended by a competent electrician.

9.2. The **Shaking Water Bath's** heating element is protected against destruction by running dry by a **low water cut-off** (thermostatical over-temperature cut-out). In case of water shortage, power is cut. A nearly unvarying, realistic actual temperature is shown on the display of the temperatur regulator, and the pilot lamps on the right side of the display and in the main switch will glow. Before restarting the unit, please fill bath with water - as described in article 4 - and unlock low water cut-off. The unlocking device is situated at the back of the unit, below a black cap nut. Within the screwing, you will see a white plastic pin, which has to be gently pressed (eg with a pen) until a clicking sound can be heard.

10. Shaking Device

1083 / 1086 Schüttelwasserbad Shaking Water Bath

The shaking device is operated by the switch in the starter panel. The shaking frequency (reciprocating) can be increased by turning the relevant knob clockwise. It can be set in a range of 10 to 250 min⁻¹.

Shaking Water Bath 1086 disposes of an LED display for the actual shaking frequency.

11. Servicing and Maintenance

GFL Shaking Water Baths are produced with first class materials and are made to withstand even rough service conditions. Nevertheless, the units should only be subjected to rough conditions within sensible limits.

Please make sure that no fluids come into contact with cable connections or the electrical part of the inside unit! Always withdraw plug before opening or cleaning the unit due to the danger of electrical shocks.

For cleaning purposes, the shaking device can be taken out of the bath. For this purpose, the device has to be pushed to the max right and unscrewed via the knurled screw (connection to motor, Pos. 4 in the exploded view). The device can now be pushed to the far left and carefully lifted out of the bath.

The inside bath can be descaled with suitable descaling agents (eg "tin-be" by Messrs. Dr. Otto Hartmann, 71665 Vaihingen/Enz, Germany).

**Do not use any hydrochlorid acid product!
These products damage the unit's interior and screwings!**

The polished exterior of the unit can easily be kept and restored with usual stainless steel polishing solutions (eg "Helios" by Messrs. Henkel Hygiene GmbH, 40589 Düsseldorf, Germany). The off-white powder coated parts can be cleaned with mild detergents if necessary.

It is advisable to regularly exchange the water within the bath. The water outlet tap is situated at the back of the unit.

The Shaking Water Bath was set and calibrated at a temperature of 50°C. For re-calibration during servicing, a calibration manual can be obtained from GFL on request. Please enquire by fax and mention type and serial number of the unit in question.

Servicings, repairs or modifications must be carried out according to the commonly recognised Technical Rules and Regulations by competent electricians only.

Only original spare parts must be used. Always demand a detailed confirmation of the carried out tasks by the person in charge (company, date, signature).

1083 / 1086 Schüttelwasserbad Shaking Water Bath

12. Technical data

Overall dimensions (wxdxh)	715x520x330 mm
Inside dimensions (wxdxh)	450x300x160 mm
Working height (incl. lid height 30 mm)	190 mm
Max. water level above shaking device / shaking tray	105 mm / 90 mm
Temperature range	
without level regulator	approx. 5°C above ambient to 99.9°C
with level regulator	approx. 3°C above tap water temperature to 99.9°C
with cooling coil (1086 only)	+10°C to 99.9°C
Temperature regulation	electronic, PI type
Temperature constancy	+/- 0.1°C (temporal)
Temperature setting and display	digital LED (0.1°C steps)
Over temperature protection	
electronic	4°C above set temperature
electro-mechanical	> 130°C via low water cut-off
Shaking motion	reciprocating
Shaking frequency	10-250 min ⁻¹
Shaking frequency display	digital LED (1 min ⁻¹ steps)
Stroke length	22 mm
Voltage	230V +/- 10%, 50....60 Hz, 1.5 kW
Power connection	Shock proof plug
Protection type / class	I / IP20
Surrounding conditions	Only inside buildings (not in explosion endangered surroundings)
Height	up to 2000 m MSL
Ambient temperature	+10°C to +40°C
Humidity	max. 80% rel. humidity, max. 31°C decreasing to 50% rel. humidity at 40°C
Weight	approx. 30 kg

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

When ordering spare parts, please always state type and serial number of the unit in question. Spare parts marked with * only refer to Shaking Water Bath 1086.

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Shaking tray made of stainless steel, perforated, to accommodate clamps for Erlenmeyer flasks of 25 ml to 500 ml. Equipped with two handles above water level for easy mounting and removal.

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Test Tube Racks made of stainless steel, equipped with two handles above water level for easy mounting and removal.

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

Flask clamps for shaking trays for Erlenmeyer Flasks, made of stainless steel, to be screwed onto Shaking Tray 3960, complete with screwing material. * = max. number of clamps per tray.

Bestell-Nr.3983 für 25 ml Kolben (52*)	Order no: 3983 for 25 ml Erlenmeyer flasks (52*)
Bestell-Nr.3984 für 50 ml Kolben (33*)	Order no: 3984 for 50 ml Erlenmeyer flasks (33*)
Bestell-Nr.3985 für 100 ml Kolben (22*)	Order no: 3985 for 100ml Erlenmeyer flasks (22*)
Bestell-Nr.3986 für 200 ml Kolben (15*)	Order no: 3986 for 200ml Erlenmeyer flasks (15*)
Bestell-Nr.3987 für 300 ml Kolben (13*)	Order no: 3987 for 300ml Erlenmeyer flasks (13*)
Bestell-Nr.3988 für 500 ml Kolben (10*)	Order no: 3988 for 500ml Erlenmeyer flasks (10*)

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Water level regulator for keeping water level constant and for cooling the Shaking Water Bath (also refer to item 5 Special Accessory Level Regulator). Adjustable, to be fixed to the back of the unit.

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