**Wasserbad  
Water Bath**

**1002 - 1013**

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**Operating  
Instructions**

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Before installation, please check whether contents of package are in good order and complete.

Should you note any damages or have any reasons for complaint, please contact your supplier or directly:

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## Contents –

1. Before installation	page	10
2. Location of the Water Bath	page	10
3. Voltage	page	10
4. Installation	page	10
5. Special Accessory Level Regulator	page	10
6. Starting Operation	page	11
7. Temperature Regulation	page	11
8. Error Display	page	12
9. Water Circulation (models 1012 and 1013 only)	page	12
10. Servicing and Maintenance	page	12
11. Technical Data	page	14
12. Circuit diagram	page	15
13. List of spare parts	page	16
14. Accessories	page	18

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## **1. Before installation**

The information given in the present manual must **by all means** be carefully read and observed. Only then can a perfect functioning of the Water Bath be guaranteed.

**Vital information within the manual are emphasized in bold letters.  
Safety precautions are additionally marked with the following symbols**

**Warning of hot surfaces**

**Warning of dangerous  
electrical voltage**

A free of charge guarantee repair cannot be granted for defects due to improper installation or handling.

## **2. Location of the Water Bath**

Place on solid, even and level surfaces only. **Do not use outside buildings.** The unit is **not** suitable for use in potentially explosive surroundings, eg during anaesthesia with inflammable gas or steam types!

## **3. Voltage**

Main switch must be in **OFF** position.

**Mains voltage and voltage stated on the name plate at the back of the unit must be identical.**

If this applies, unit can be connected.

**The Water Bath must only be connected to a correctly installed shock proof power socket.**

## **4. Installation**

**The Water Bath must only be operated with tap water. Other media, such as eg oil, acid or any agents not approved by GFL can lead to damages of the bath interior, the screw connections or the heating element.**

Before operation, the bath interior must be filled with water. Do make sure that the water outlet tap at the back of the unit is closed (handle must have a 90° angle from tap position). Open the lid and fill the bath. The water level should be kept between the "**min**" (50 mm above the perforated tray and "**max**" (50 mm below the upper bath rim) markings.

## **5. Special Accessory Level Regulator**

With the adjustable level regulator at the back of the unit, occurring water loss (through evaporation) is compensated and thus water level is kept constant. The level regulator also allows operation of the bath below ambient, minimum of approx. 3°C above tap water temperature. Before operation, the level regulator must be connected to water tap.


## 1002 - 1013 Wasserbäder Water Baths

---

The upper hose connection of the water level regulator is the water supply. Connect this to the tap a laboratory hose (inside diameter max. 9 mm). The lower hose connection is for water outlet and must be connected to a lower placed drain. Make sure that the water flows off without being dammed up.

**All hose connections must be secured with hose clamps.**

The required water level can be adjusted with the outlet tube of the level regulator. To do so, loosen the outlet tube screwing with a spanner GW 27. The water level can now be set by pulling out or pushing in the tube. Tighten screwing again. Now carefully open water tap and fill Water Bath with water as described in article 4.

### 6. Starting Operation

Switch on main switch. The green pilot lamp will glow.

### 7. Temperature Regulation

After switching on the **Water Bath**, the display shows the actual inside basin temperature. By pressing key "°C", the display will show the last set and registered temperature that will now be used again automatically.

The temperature range is between approx. 5°C above ambient to 99.9°C, temperature constancy (temporal) approx. +/- 0.1°C. The set temperature can be changed with keys "°C", "+" and/or "-". By pressing key "°C", the display switches from actual to set temperature. By pressing key "°C" and either "+" and/or "-" simultaneously (two finger operation), the required temperature can be set in 0.1°C steps. The counting speed increases after a short time. The new set temperature is saved as soon as key "°C" is released. The display will now show the actual temperature again.

The yellow pilot lamp next to the temperature display shows that heating is in operation.

**CAUTION:**  
**Hot surfaces when temperatures  
rise above 60°C .  
Do not touch!  
Danger of burns!**

## **8. Error Display**

8.1. The temperature regulator disposes of electronic monitoring with error screening and concurrent over-temperature cut-out.

The **Temperature Sensor** of the regulator is constantly being screened for short circuits. In case of fault, display will show: "**E1L**" for a short circuit in the Temperature Sensor or "**E1H**" for interruptions. The microprocessor controlled temperature regulator will switch off the heating circuit.

After elimination of error, unit can be restarted.

The unit also disposes of an **electronic over-temperature cut-out** for the protection of test substances, dependent on the set temperature. This device will switch off heating when the set temperature is exceeded by 4°C. The electronic monitoring will induce the information "**E2H**" to be shown on the display.

To restart heating, the unit must be switched off and then on again at the main switch. Should the information "E2H" remain, error must be attended by a competent electrician.

8.2. The **Water Bath's** heating element is protected against destruction by running dry by a **low water cut-off** (thermostatical over-temperature cut-out). In case of water shortage, power is cut. A nearly unvarying, realistic actual temperature is shown on the display of the temperature regulator, and the pilot lamps on the right side of the display and in the main switch will glow. Before restarting the unit, please fill bath with water - as described in article 4 - and unlock low water cut-off. The unlocking device is situated at the back of the unit, below a black cap nut. Within the screw connection, you will see a white plastic pin, which has to be gently pressed (eg with a pen) until a clicking noise can be heard.

## **9. Water Circulation (models 1012 and 1013 only)**

The integrated water circulation system ensures optimum temperature distribution within the whole basin. It is operated by an electric motor with stirring magnet flanged to the bath basin. Its torque is transmitted to a teflon coated stirrer, which draws water from the middle of the basin and then pushes the water back to all sides. The circulation motor is activated via the main switch.

## **10. Servicing and Maintenance**

GFL Water Baths are produced with first class materials and are made to withstand even rough service conditions. Nevertheless, the units should only be subjected to rough conditions within sensible limits.

**For technical support concerning GFL Water Baths,  
please do not hesitate to contact our technical staff.**

Before returning a defective unit, please contact your dealer or us.

**Please do not return the unit before consulting us.**

If we agree to the unit being returned, arrange for careful packing and send the unit to

**GFL Gesellschaft für Labortechnik mbH**  
**Schulze-Delitzsch-Strasse 4**  
**D-30938 Burgwedel**  
**Federal Republic of Germany**

**Please make sure that no fluids come into contact with cable connections or the electrical part of the inside unit!  
Always pull the plug before opening or cleaning the unit due to the danger of electrical shocks.**

For cleaning purposes, the perforated tray can be taken out of the bath.  
The inside bath can be descaled with suitable descaling agents (eg "tin-be" by Messrs. Dr. Otto Hartmann, 71665 Vaihingen/Enz, Germany).

**Do not use any hydrochlorid acid product!  
These products damage the unit's interior and screw connections!**

The polished exterior of the unit can easily be kept and restored with usual stainless steel polishing solutions (eg "Helios" by Messrs. Henkel Hygiene GmbH, 40589 Düsseldorf, Germany). The off-white powder coated parts can be cleaned with mild detergents if necessary.

It is advisable to regularly exchange the water within the bath. The water outlet tap is situated at the back of the unit.

The GFL Water Bath was set and calibrated at a temperature of 50°C. For re-calibration during servicing, a calibration manual can be obtained from GFL on request. Please enquire by fax and mention type and serial number of the unit in question.

Servicings, repairs or modifications must be carried out according to the commonly recognised Technical Rules and Regulations by competent electricians only.

Only original spare parts must be used. Always demand a detailed confirmation of the carried out tasks by the person in charge (company, date, signature).

# 1002 - 1013 Wasserbäder Water Baths

---


## 11. Technical Data

Model/overall dimensions (wxdxh in mm)/ Inside dimensions (wxdxh in mm)	1002/325x395x255 / 245x200x145 1003/485x440x255 / 400x245x145 1004/685x440x255 / 600x245x145 1005/495x490x445 / 410x296x326 1008/485x440x325 / 400x245x205 1012/325x395x325 / 245x200x145 1013/485x440x325 / 400x245x145
Model/working height (incl. lid 30 mm)/ Max. water level above perforated tray/ volume (in litres)	1002 / 175 mm / 95 mm / 7 l 1003 / 175 mm / 95 mm / 14 l 1004 / 175 mm / 95 mm / 21 l 1005 / 345 mm / 265 mm / 40 l 1008 / 235 mm / 160 mm / 20 l 1012 / 175 mm / 95 mm / 7 l 1013 / 175 mm / 95 mm / 14 l
Temperature range without level regulator 99.9°C with level regulator	approx. 5°C above ambient to  approx. 3°C above tap water temperature to 99.9°C
Temperature regulation Temperature constancy Temperature setting and display Over temperature protection electronic electro-mechanical	electronic, PI type +/- 0.1°C (temporal) digital LED (0.1°C steps)  4°C above set temperature > 130°C via low water cut-off
Voltage Model / power	120V or 230V +/- 10%, 50....60 Hz 1002 / 1.0 kW 1003 / 1.5 kW 1004/ 1.5 kW 1005 / 1.5 kW 1008 / 1.5 kW 1012 / 1.0 kW 1013 / 1.5 kW
Power connection Protection type / class Surrounding conditions	Shock proof plug I / IP20 Only inside buildings (not in potentially explosive surroundings) up to 2000 m MSL +10°C to +40°C max. 80% rel. humidity max. 31°C, decreasing to 50% rel. humidity at 40°C
Height Ambient temperature Humidity	
Model / weight	1002 / 8.0 kg 1003 / 12.0 kg 1004 / 15.0 kg 1005 / 21.0 kg 1008 / 13.0 kg 1012 / 8.0 kg 1013 / 12.0 kg

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
When ordering spare parts, please **always state model and serial numbers** of the unit in question.

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**Water Level Regulator** for keeping water level constant and for cooling the Water Bath (also refer to item. 5 Special Accessory Level Regulator). Adjustable, to be fixed to the back of the unit.

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**Stainless Steel Racks** for test tubes and bottles

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Bestell-Nr. 1912

6 Flaschen (6x200 ml)

ProAquaTop prevents the  
formation of algae, bacteria  
and mould in the Water Bath.  
With colour indicator to  
control the concentration of  
the solution of 0.1%.

Order No. 1910

1 Bottle (200 ml)

Order No. 1911

3 Bottles (3x200 ml)

Order No. 1912

6 Bottles (6x200 ml)

**Gesellschaft für Labortechnik mbH**  
Postfach / P.O. Box 1152 · D-30927 Burgwedel  
Schulze-Delitzsch-Straße 4 · D-30938 Burgwedel  
Telefon +49 5139 9958 0 · Telefax +49 5139 9958 21  
E-Mail: [info@GFL.de](mailto:info@GFL.de)  
Internet: <http://www.GFL.de>