*illumipro-10*TM

Operator's Manual

 **Meridian**
Bioscience, Inc.
Inspired Science. Trusted Solutions.®

Table of Contents

| | |
|---|--------------|
| Instrument Intended Use and Function | 3 |
| Installation Procedures and Requirements | 3 |
| Principles of Operation | 5 |
| Performance Characteristics and Specifications | 5 |
| Operating Instructions | 7 |
| Calibration Procedures | 13 |
| Operational Precautions and Limitations | 15 |
| Service and Maintenance | 16 |
| <i>illumipro-10</i> Set-Up and Installation Checklist | Appendix I |
| <i>illumipro-10</i> Decontamination Record | Appendix II |
| <i>illumipro-10</i> Error Codes and Troubleshooting | Appendix III |
| <i>illumipro-10</i> Barcode Scanner Information | Appendix IV |
| <i>illumipro-10</i> Compliance Testing Summary | Appendix V |

Instrument Intended Use and Function

The *illumipro-10* is an automated isothermal amplification and detection system for use with Meridian Bioscience, Inc. *illumigene* Loop-Mediated Amplification products.

The *illumipro-10* is intended to be used by trained laboratory professional users in laboratory settings.

Installation Procedures and Requirements

The *illumipro-10* and its accessories are securely packaged to prevent damage during shipment to the end user. The *illumipro-10* shipping container and packaging should be inspected for damage prior to installation. Damaged instrumentation should not be installed as this could create a hazard to the end user. Report any damage to Meridian's Technical Support staff at 1-800-343-3859.

illumipro-10 Package Contents

- ① *illumipro-10* Automated Isothermal Amplification and Detection System
- ② *illumipro-10* Power Supply and Power Cord
- ③ *illumipro-10* Verification Standards
- ④ *illumipro-10* Operator's Manual Binder
- ⑤ *illumipro-10* USB Cable

External Thermal Printer Package Contents (*Catalog 610173; shipped separately*):

- ① Printer
- ② Printer Cable
- ③ Printer Power Supply (12VDC Linear Adaptor)
- ④ Thermal Roll Paper (1 Roll)

illumipro-10 Optional Accessories

- ① External Keyboard (*Catalog 610174, shipped separately.*)

illumipro-10 Operator's Manual (*SN11007; shipped separately.*)

illumipro-10 Installation

Installation of the *illumipro-10* and its accessories can be performed after contents have been inspected and the requirements described in this Operator's Manual have been reviewed. Appendix I includes a general Checklist for Set-up and Installation of the *illumipro-10*.

The *illumipro-10* should be placed on a sturdy, level surface. Set-up the instrument for use as described in the table below.

| Component | Symbol | Installation Instruction |
|--------------------------|---|--|
| Printer |  | <p>Connect the external printer to the <i>illumipro-10</i> using the supplied cabling. The printer connection is located at the back of the instrument and is identified by the symbol shown.</p> <p>Secure the printer connection and connect the power supply cable provided with the printer packaging to the printer.</p> |
| Power Supply | 12V ---- @4.5A | <p>Connect the power supply cable provided with the <i>illumipro-10</i> packaging to the unit. The power supply connection is located at the back of the instrument. The connection port for the power supply is identified by the symbol shown. Connect the power supply cable to the power supply box.</p> <p>Plug the pronged ends of the <i>illumipro-10</i> and printer power cords to an appropriate power receptacle.</p> |
| External Keyboard |  | <p>Connect the optional, external keyboard to the unit using the supplied cabling. The External Keyboard connection is located at the back of the instrument and is identified by the symbol shown.</p> <p>NOTE: The Keyboard must be installed with the <i>illumipro-10</i> powered off.</p> |

illumipro-10 Set-up

Instrument set-up is completed under the SYSTEM Menu. The user will be able to set the time format, date format and preferred language.

illumipro-10 Performance Verification

Performance verification must be performed after installation and prior to use. Optical verification is completed according to the instructions provided under the SERVICE Menu.

Principles of Operation

The *illumipro-10* is an automated isothermal amplification and detection system for target nucleic acid sequences found in human specimens. The instrument is used in conjunction with Meridian Bioscience, Inc's *illumigene* LOOP-Mediated Amplification in vitro diagnostic products.

The *illumipro-10* is a menu driven laboratory instrument with two independent sample processing blocks, identified as Block A and Block B. Sample heating and optical detection is carried out for up to five two-chambered *illumigene* devices per Block. Each two-chambered *illumigene* device contains a Sample Chamber and a Control Chamber. Amplification of target DNA occurs during the heat cycle and results in the formation of precipitate detected by the *illumipro-10* optics system. The precipitate generated by the presence of amplified target DNA leads to a turbid Sample/Control reaction solution which is then measured by absorbance. The *illumipro-10* uses the change in turbidity of each Sample/Control reaction solution to report assay results as INVALID, POSITIVE, or NEGATIVE.

The *illumipro-10* operates in four basic modes: ASSAY, RESULTS, SERVICE, and SYSTEM. Assay Selection and Sample Amplification occurs in the ASSAY mode; Test Results are managed in the RESULTS mode; basic instrument set-up is performed in the SYSTEM mode; and optical performance verification is completed in the SERVICE mode.

Performance Characteristics and Specifications

Performance Characteristics

The *illumipro-10* is an automated isothermal amplification and detection system for use with Meridian Bioscience, Inc. *illumigene* Loop-Mediated Amplification products. The *illumipro-10* was designed with a simple User Interface which includes a keypad, liquid crystal display (LCD), barcode scanner, printer and optional external keyboard. Menu-driven operating instructions are shown on the display and the user inputs commands to the instrument by making selections through the keypad.

Isothermal amplification is carried out by two independently controlled heat blocks capable of operating between 55 C and 65 C, and within 1 C of the temperature set point. The isothermal amplification temperature set point is dictated by the *illumigene* Assay selected. Isothermal amplification time is monitored by the *illumipro-10*'s internal timer. When in operation, the *illumipro-10* will display 'Test in Progress' as well as block temperature and incubation time remaining.

DNA amplification detection is completed by the *illumipro-10* optics system. Each *illumipro-10* block contains laser diodes illuminating at 650 ± 20 nm and corresponding detectors monitoring light transmission across each *illumipro-10* well. The *illumipro-10* performs an initial check of the optics system prior to run start. Observed failures in the optics system disable the instrument block until the failure can be resolved. After successful optics verification, the *illumipro-10* checks for the presence of an *illumigene* test device in each well. The *illumipro-10* automatically enters a Sample ID of 'EMPTY WELL' for all wells that a test device cannot be detected in. The *illumipro-10* measures the absorbance of each *illumigene* test device at the beginning and end of the isothermal amplification incubation. Sample results are reported as INVALID, POSITIVE, or NEGATIVE based on the observed change in absorbance.

Performance features of the *illumipro-10* include the following:

- Self-diagnostics are performed at power-on and at Run start.
- Self-diagnostic for Test Device Placement and Optics Path Verification.
- Lid Closure Detection and Lid Lock Feature.
- Thermal Heat Block with Optics Detection at each well with Visible Light Emission/Absorbance (650 ± 20 nm).
- Interval Timer with Real-Time Assay Processing Time Display.
- Simple LCD User Interface for Instrument Set-up, Program Selection and Sample Identification Input.
- Barcode Scanner for Sample Identification entry.
- External Keyboard available for Sample Identification entry.
- Attached Printer for results reporting.
- Control of amplification product carry-over through use of the self-contained *illumigene* Test Device.

illumipro-10 Specifications

- **Electrical**

| | |
|-----------------------------|--------------------------|
| Voltage Range: | 120 V AC |
| Operating Range, Supply: | 100 – 240 V AC, 50/60 Hz |
| Voltage and Current Rating: | 12 V DC, 4.5 Amp |
- **Physical**

| | |
|-------------|--------------------------|
| Dimensions: | 21 cm x 29.2 cm x 9.5 cm |
| Weight: | 2.95 \pm 0.05 kg |
- **Environmental**

| | |
|-------------------------------|--------------------------|
| Operating Temperature: | 15 – 30 C |
| Storage Temperature: | 10 – 40 C |
| Relative Humidity, Operating: | 10 – 90%, non-condensing |
| Relative Humidity, Storage: | 10 – 95% |

Printer Specifications

- **Electrical**

| | |
|------------------------|------------------------------|
| Supply Voltage (SMPS): | Input voltage (S, P) 12 V DC |
|------------------------|------------------------------|
- **Environmental**

| | |
|-------------------------------|------------|
| Operating Temperature: | 0 – 40 C |
| Storage Temperature: | -10 – 50 C |
| Relative Humidity, Operating: | 30 – 80% |
| Relative Humidity, Storage: | 10 – 90% |

Operating Instructions

The *illumipro-10* operates in four basic modes: ASSAY, RESULTS, SERVICE, and SYSTEM. Assay Selection and Sample amplification occurs in the ASSAY mode; Test Results are managed in the 'RESULTS mode; and basic instrument set-up is performed in the SYSTEM mode. The SERVICE mode is reserved for trained service professionals and is not accessible by the Laboratory User. General information regarding each mode of operation is provided in this section.

Keypad

The functions of the *illumipro-10* are navigated through the keypad. The Keypad provides a simple user interface and allows for basic menu navigation, input of alphanumeric characters for sample identification, and Assay RUN initiation. Keypad functions will be referred to throughout this manual; keypad layout and symbols used are defined below.


| Keypad Button | Character / Function |
|---------------|-----------------------------|
| 1 | <SPACE>, 1 |
| 2 | A, B, C, 2 |
| 3 | D, E, F, 3 |
| 4 | G, H, I, 4 |
| 5 | J, K, L, 5 |
| 6 | M, N, O, 6 |
| 7 | P, Q, R, S, 7 |
| 8 | T, U, V, 8 |
| 9 | W, X, Y, Z, 9 |
| 0 | 0 |
| - | - |
| ▲ | Scroll Up |
| ▼ | Scroll Down |
| ◀ | Scroll Back, Backspace |
| ↵ | Enter |
| RUN | Run Selected Assay Protocol |

NOTE: For keypad buttons with multiple characters, scroll by pressing the keypad button multiple times.

Main Menu

The Main Menu screen allows the user to view Time Date, Assay Status and Block Temperature. The Main Menu is used to access the Results, Service and System Modes.


Conventions Used:

| Graphic | Description |
|---|---|
| *** | BLOCK STATUS: Warming up; Block not at temperature |
| --:-- | BLOCK STATUS: Idle; Block at Temperature |
| !!! | Warning, Instrument Requires Attention |
| ... | Assay Menu Expansion Indicator |
|  | Printer: Connect Printer to unit at this port. (Back of instrument) |
|  | USB Port: Connect unit to external computer at this port. (Back of instrument) |
| PS/2  | External Keyboard: Connect unit to external keyboard (optional) at this port. (Back of instrument) |
| 12V ----- @4.5A | illumipro-10 Power Supply: Connect Power Supply to unit at this port. (Back of instrument) |


ASSAY MODE

The ASSAY MODE allows the user to access and run programs on the *illumipro-10*. The user selects the Block to be used from the Main Menu and then selects the Assay to be run. Each Block can be selected and run independently. Upon selection of the Block to be used, the Assay Mode Menu appears. The User selects the assay to be run and follows the instructions displayed on the *illumipro-10* display.

The *illumipro-10* display indicates that a test is in process and prompts the user to enter Sample Identification information or to abort the test. Sample Identification information can be entered directly using the Keypad, may be scanned using the *illumipro-10* barcode scanner or may be entered using the external keyboard. Additional information regarding the *illumipro-10* barcode scanner is provided in Appendix IV of this manual.


The User enters Sample Identification Information by following the instructions shown on the *illumipro-10* display and initiates the Assay Run.

NOTE: Block temperature will appear in the 'TEMP' field. The 'Block Status' field will display a timer indicated time until assay completion.

The *illumipro-10* will complete the selected assay program and the Block Status on the MAIN MENU will be shown as '00:00'. Results for the selected assay can be viewed by selecting the Block displaying completed status and following the instructions shown on the *illumipro-10* display. Results can be printed manually by selecting the print option at the end of the results display or by enabling the Auto-Print feature during System Set-up.


Upon completion of the assay run, the user must remove the *illumigene* devices from the unit and discard appropriately.


NOTE: *Care must be taken to avoid contamination of the equipment and workspace by target and/or amplified nucleic acids. DO NOT open illumigene devices after assay completion.*

RESULTS MODE

The RESULTS MODE allows the user to view stored results, delete stored results and enable the auto-printer. The *illumipro-10* will store up to 1000 individual test results or 200 Batches. The instrument displays a warning to the user when result storage is approaching maximum capacity. RESULTS MODE is accessed by following the instructions shown on the *illumipro-10* display.


The User can view stored results by selecting the 'View Results' option. The 'View Results' menu allows the user to access stored data by Date and Block. Selected Batch Data is displayed and can be printed manually by selecting the print option at the end of the results display.

SERVICE MODE


The SERVICE MODE allows the user to perform optics system verification, view assay parameters, print system check information, and configure the printer. **NOTE:** The SERVICE MODE menu includes a 'Diagnostics' Option which can be accessed only by trained service personnel.

OPTICS SYSTEM VERIFICATION is required to ensure proper function of the *illumipro-10*. Instructions for completion of the Optics System Verification are provided within the **CALIBRATION PROCEDURE** Section of this manual.

SYSTEM MODE

The SYSTEM MODE allows the user to set and format time and date, configure the printer and barcode scanner, set user language, and view system information. The user may set time, date, or language (English, Italian, French, Spanish and German) based on local requirements and/or preference. The SYSTEM MODE allows the user to enable or disable automatic printing and the barcode scanner. System configuration is completed by following the instructions shown on the *illumipro-10* display.

NOTE: *Time and Date cannot be modified while a test is in process. Language cannot be modified while a test is in process.*



Calibration Procedures

POWER-ON SELF TEST

The *illumipro-10* automatically performs an internal **Power-On Self Test (POST)** each time the instrument is powered-on. The POST Test confirms that software and hardware components of the system are performing as expected. Successful completion of the POST Test will be followed by an audible tone. POST Test failure will be indicated by an Error Code. Additional information regarding Error Codes is located in Appendix III of this manual.

OPTICS SYSTEM VERIFICATION

Calibration of the *illumipro-10* is not required. Verification of the OPTICS SYSTEM for each Block must be performed monthly to ensure proper function. Verification is completed using the Red Verification Standard included with the *illumipro-10*. Stepwise instructions for Optics System Verification are provided on the *illumipro-10* display, utilizing the SERVICE MODE menu.

Optic System Verification includes two stages. The first stage ensures that the optical path is clear and free of obstruction. The second stage confirms that the optics system is transmitting and detecting emitted light particles properly. Transmission and detection verification requires the use of the *illumipro-10* Red Verification Standard. The *illumipro-10* steps through the Verification protocol and prompts the user for action as necessary.

The Red Verification Standard must be seated firmly in the Heat Block wells. The Verification Standard should be oriented with the asset number tag at the Well 1 position.

Upon completion of the Verification Protocol, the *illumipro-10* will display verification test results as either 'Pass' or 'Fail'. In the event that Verification testing does not give acceptable results, the user should verify the optical path is free of obstructions, the standard is free of visible defects and repeat the verification testing. If repeat testing does not give passing verification results, contact Meridian's Technical Support Staff for further assistance.

NOTE: Each Block of the *illumipro-10* functions independently. For example, a failed verification test for Block A will not prevent use of Block B when Block B verification testing passes.

WARNING: *Do NOT leave verification standards in the illumipro-10. Verification Standards will get hot and could burn the user. Exposure to temperature may impact the performance of the Verification Standards.*

Operational Precautions and Limitations

WARNINGS:

| | |
|---|---|
|  | <p>CAUTION: Risk of Danger. The <i>illumipro-10</i> is an electromechanical device that can cause physical shock or injury to the operator if not used in accordance with the procedures described in this manual.</p> |
|  | <p>LASER RADIATION: Avoid Exposure to Beam. The <i>illumipro-10</i> contains a class 3R laser product. The laser will not function when the lid is in the open position, however, care should be taken in the handling and use of this instrument.</p> |
|  | <p>HOT SURFACE: Keep hands Away from Hot Surfaces. The <i>illumipro-10</i> contains a heat block producing temperatures between 55 – 65 C during operation. Care should be taken to avoid direct contact with the heat block.</p> |
|  | <p>CAUTION: Laser Radiation. The <i>illumipro-10</i> contains a laser product. Service of the unit should be performed only by qualified personnel as optical exposure to the laser beam may cause injury.</p> |
| <p>IPX-0</p> | <p>CAUTION: Protect from water. The <i>illumipro-10</i> does not protect from the egress of water. Do not expose to water or submerge the instrument.</p> |

PRECAUTIONS:

- The *illumipro-10* is an automated instrument that utilizes isothermal Loop-Mediated Amplification Technology. Care must be taken to avoid contamination of the equipment and workspace by target and/or amplified nucleic acids. Only qualified personnel should perform molecular testing.
- The *illumipro-10* is used in conjunction with Meridian Bioscience, Inc's *illumigene* Loop-Mediated Amplification in vitro diagnostic products. Samples processed in the *illumipro-10* should be handled in accordance with the instructions provided in specific *illumigene* product instructions for use.
- Selected Language cannot be changed while a run is in process.
- Time and Date cannot be changed while a run is in process
- Optical System Verifiers must be stored in the case provided. The optical verifiers must be protected from light and damage. Optical system verifiers should not be stored in the *illumipro-10*.
- Do not hot swab the optional *illumipro-10* keyboard. CAPS Lock will come on at power-up to indicate keyboard is online.
- The optional Keyboard must be plugged in while the *illumipro-10* is powered-off.
- When not in use, the *illumipro-10* should be stored with the lid closed.

Service and Maintenance

Servicing of the *illumipro-10* should be performed by qualified professionals only. Contact Meridian Bioscience Technical Support at (800) 343-3858 for technical support or to make arrangements for service. DO NOT ATTEMPT TO SERVICE THE *illumipro-10*.

Surface Cleaning

Cleaning of the exterior surfaces of the *illumipro-10* and the immediate work area should be performed as necessary, no less than daily when in use. Allow instrument to cool and wipe surfaces with a lint-free cloth moistened with DNase/RNase cleaning solution or a 10% Bleach Solution.

WARNING: Surface cleaning should be performed only when the instrument is powered-off AND disconnected from the power source. DO NOT use saturated clothes for cleaning.

Heat Block Cleaning

Cleaning of the instrument heat block should be performed by qualified personnel only. Cleaning of the heat block should be performed only when contamination of the heat block is suspected. Contamination of the heat block may be DNA Sourced, or non-DNA Sourced. The cleaning protocol followed should be based on the source of the contamination as shown below:

- **DNA Contamination Cleaning Protocol**
 1. Gently wipe heat block chamber with a dry Foam Swab.
 2. Gently wipe heat block chamber with a Foam Swab moistened with 10% Bleach Solution.
 3. Gently wipe heat block chamber with a dry Foam Swab.
 4. Gently wipe heat block chamber with a Foam Swab moistened with 70% Alcohol.
 5. Gently wipe heat block chamber with a dry Foam Swab.

DO NOT use saturated Foam Swabs for cleaning.

- **non-DNA Contamination Cleaning Protocol**
 1. Gently wipe heat block chamber with a dry Foam Swab.
 2. Gently wipe heat block chamber with a Foam Swab moistened with 70% Alcohol.
 3. Gently wipe heat block chamber with a dry Foam Swab.

DO NOT use saturated foam swabs for cleaning.

WARNING: Heat block cleaning should be performed only when the instrument is powered-off AND disconnected from the power source. DO NOT use saturated swabs for cleaning.

WARNING: Do not attempt to clean the *illumipro-10* using compressed air.

ALWAYS perform Optical Verification Testing after Heat Block Cleaning.

Appendix I

illumipro-10 Set-Up and Installation Checklist

| Facility Information | | | | |
|----------------------------|---|------------------------------------|--------------------------------|---------------------------------|
| Facility Name: | | | | |
| Facility Address: | | (Street Address) | | |
| | | (City, State) | | |
| Instrument Information | | | | |
| Instrument Serial Number: | | | | |
| Installation Location: | | (Room, Floor, Building, et cetera) | | |
| Installation Performed By: | | | Installation Date: | |
| Set-Up and Installation | | | | |
| 1. | The <i>illumipro-10</i> shipping container and packaging have been inspected for damage; no damage found. Installation of damaged could create a hazard to the end user. | Yes <input type="checkbox"/> | No <input type="checkbox"/> | N/A <input type="checkbox"/> |
| 2. | The <i>illumipro-10</i> Operator's Manual has been received and reviewed. Operational Precautions and Limitations have been reviewed and are understood. | Yes <input type="checkbox"/> | No <input type="checkbox"/> | N/A <input type="checkbox"/> |
| 3. | The Installation Location meets Electrical and Environmental Specifications described in the <i>illumipro-10</i> Operator's Manual. | Yes <input type="checkbox"/> | No <input type="checkbox"/> | N/A <input type="checkbox"/> |
| 4. | Installation of the <i>illumipro-10</i> and its accessories has been completed as described by the Operator's Manual. | Yes <input type="checkbox"/> | No <input type="checkbox"/> | N/A <input type="checkbox"/> |
| 5. | The <i>illumipro-10</i> has been powered-on and the Power On Self-Test was successful. | Yes <input type="checkbox"/> | No <input type="checkbox"/> | N/A <input type="checkbox"/> |
| 6. | Set-up of the <i>illumipro-10</i> including Time Format, Date Format, and Preferred Language has been completed as described in the Operator's Manual. | Yes <input type="checkbox"/> | No <input type="checkbox"/> | N/A <input type="checkbox"/> |
| 7. | <i>illumipro-10</i> Performance Verification was completed as described in the Operator's Manual. Acceptable results were obtained for Block A and Block B. | Yes <input type="checkbox"/> | No <input type="checkbox"/> | N/A <input type="checkbox"/> |
| 8. | A copy of the printout for the System Check and Optics System Verification is attached to this record. | Yes <input type="checkbox"/> | No <input type="checkbox"/> | N/A <input type="checkbox"/> |
| Installation Complete | | | | |
| (Signature/Date) | | | | |

Appendix II

illumipro-10 Decontamination Record

Decontamination of laboratory equipment must be completed prior to return of the instrument for servicing. Decontamination must be completed according to the instructions provided. Please complete the information below and attach the completed record to the instrument prior to return.

Instruments returned for service without evidence of decontamination will not be serviced. If the instrument is submitted for service without decontamination documentation, it will be shipped back to the user at the user's expense.

1. Wipe exterior surfaces of the *illumipro-10* with a lint-free towel moistened with a 10% Bleach-based Disinfecting Solution. Minimum contact exposure time for kill of Bloodborne Pathogens (Hepatitis-A, HIV-1, MRSA, SARS, et cetera) is 1 minute.
2. Open *illumipro-10* Lids A and B. Wipe inside lid surface and Heat Block Surfaces with lint-free towel moistened with 10% Bleach-based Disinfecting Solution. Minimum contact exposure time is 1 minute.
3. Close Lids. Affix the completed *illumipro-10* Decontamination Record to instrument. Package the instrument in the protective foam packaging and shipping container provided.

| Instrument Information | |
|---|------------------|
| Instrument Serial Number: | |
| Facility Information | |
| Facility Name: | |
| Facility Address: | (Street Address) |
| | (City, State) |
| Contact Information: | (Name, Title) |
| | (Telephone) |
| | (e-mail) |
| Decontamination Information | |
| The <i>illumipro-10</i> has been decontaminated as described above. Exterior and Interior contact surfaces have been wiped with a 10% Bleach-based Disinfecting Solution. Minimum contact exposure times have been monitored. | |
| Performed By: | (Signature/Date) |

Appendix III

illumipro-10 Error Codes and Troubleshooting

The *illumipro-10* will display Error Codes when system failures are detected. Error Codes, descriptions and the actions recommended for correction are summarized in the table below.

| Error Code | Description | Action Required |
|------------|--|---|
| 500 | Real time clock not working. | Contact Meridian Bioscience Technical Support at (800) 343-3858 for technical support or to make arrangements for service. |
| 501 | Non-Volatile RAM (NVRAM) not working. | |
| 504 | Voltage out of range. | |
| 505 | Voltage out of range. | |
| 506 | Voltage out of range. | |
| 507 | Voltage out of range. | |
| 508 | Interior temperature is too hot. | Move unit to a cooler environment |
| 509 | Block A temperature is too hot. | Turn unit off and wait for it to cool. |
| 510 | Block B temperature is too hot. | |
| 502 | Firmware Error: Heater block is not heating properly. | Power-off instrument. |
| 503 | Firmware Error: Non-Volatile RAM (NVRAM) did not initialize properly. | Turn instrument back on and allow POST Testing to be completed. Contact Meridian Bioscience Technical Support at (800) 343-3858 for technical support if problem persists. |
| 511 | Firmware Error: Test results were not fetched from Non-Volatile RAM (NVRAM) properly. | |
| 512 | Firmware Error: Device parameters were not fetched from Non-Volatile RAM (NVRAM) properly. | |
| 513 | Firmware Error: Block status cannot be read properly. | |

Appendix IV

illumipro-10 Barcode Scanner Information

The *illumipro-10* allows for entry of Sample identification using barcodes. The barcode scanner has been placed on front of the instrument. The scanner is oriented vertically, next to the keypad.

Sample Identification information is limited to 16 characters. Barcodes containing only the following characters will be accepted:

| | | | | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|---|---|---|
| A | B | C | D | E | F | G | H | I | J | K | L | M |
| N | O | P | Q | R | S | T | U | V | W | X | Y | Z |
| 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | - | | |

Appendix V

illumipro-10 Compliance Testing Summary

The *illumipro-10* has been tested and found to be in compliance with the following requirements and Standards:

| Standard | Description |
|-----------------------|---|
| IEC 61010-1 | Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements |
| IEC 61010-2-010 | Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials |
| IEC 61010-2-081 | Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes |
| IEC 61010-2-101 | Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment |
| UL 61010-1 | UL standard for Safety Electrical Equipment for Measurement, Control, and Laboratory Use - Part 1: General Requirements |
| CSA C22.2#61010-1 | UL standard for Safety Electrical Equipment for Measurement, Control, and Laboratory Use - Part 1: General Requirements |
| CSA C22.2#61010-2-010 | Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials |
| CSA C22.2#61010-2-081 | Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use - Part 2-081: Particular Requirements for Automatic and Semi-Automatic Laboratory Equipment for Analysis and Other Purposes |
| CSA C22.2#61010-2-101 | Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory use - Part 2-101: Particular Requirements for In Vitro Diagnostic (IVD) Medical Equipment |

This page intentionally left blank.

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Anlage IV

illumipro-10 Barcode-Scanner-Informationen

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| Standard | Description |
|-----------------------|---|
| IEC 61010-1 | Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements |
| IEC 61010-2-010 | Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials |
| IEC 61010-2-081 | Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes |
| IEC 61010-2-101 | Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment |
| UL 61010-1 | UL standard for Safety Electrical Equipment for Measurement, Control, and Laboratory Use - Part 1: General Requirements |
| CSA C22.2# 61010-1 | UL standard for Safety Electrical Equipment for Measurement, Control, and Laboratory Use - Part 1: General Requirements |
| CSA C22.2#61010-2-010 | Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials |
| CSA C22.2#61010-2-081 | Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use - Part 2-081: Particular Requirements for Automatic and Semi-Automatic Laboratory Equipment for Analysis and Other Purposes |
| CSA C22.2#61010-2-101 | Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory use - Part 2-101: Particular Requirements for In Vitro Diagnostic (IVD) Medical Equipment |

SN11007

REV. 07/10

| | |
|--|---|
|  Manufactured By | <p>Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: (513) 271-3700 Orders/Customer Service: (800) 543-1980 Technical Support Center: (800) 343-3858 Information Fax:(513) 272-5432 Ordering Fax:(513) 271-0124</p> |
|  Authorized Representative | <p>Meridian Bioscience Europe Via dell' Industria, 7 20020 Villa Cortese (MI) Italy Tel.: +39 0331 433636 Fax: +39 0331 433616 e-mail: info@mdeur.com</p> |

Meridian Bioscience Europe France
Le Quadra
455, Promenade des Anglais
06299 Nice Cedex-3
France
Tel.: +33 (4) 93 18 72 10
Fax: +33 (4) 93 18 72 11
e-mail: info@meridianbioscience.fr

Meridian Bioscience Europe s.a. / n.v.
Rue de l'Industrie 7
1400 Nivelles
Belgium
Tel.: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
e-mail: info@mdeur.be

Meridian Bioscience Europe b.v.
Halderheiweg 6
5282 SN Boxtel
The Netherlands
Tel.: +31 (411) 621166
Fax: +31 (411) 624841
e-mail: meridian.info@planet.nl

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symbols, Guia de simbolos, Erläuterung der graphischen symbole)

| | | | |
|-------------|--|-------------------------|--|
| | Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis | CONTROL + | Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle |
| LOT | Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung | CONTROL - | Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle |
| IVD | In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum | SMP PREP DIL SPE | Sample Preparation Apparatus containing Sample Diluent / Apparato per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet |
| CE | This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG. | EC REP | Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft |
| REF | Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer | | Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren |
| | Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten | RoHS | Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe |
| | Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller | | Caution, consult accompanying documents / Attenzione, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten |
| | Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen | BUF RXN | Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer |
| | Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung | | ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt |
| SN | Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer | | Recycle - do not dispose of as general waste / Riciclare - non eliminare come rifiuto generico / Recycler - ne pas jeter dans une poubelle / Recycle - no desecha como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen |
| TEST | Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät | EX TUBE | Extraction tube / Provette per l'estrazione / Tubes d'Extraction / Tubo de Extracción / Röhrchen zur Probenextraktion |
| | Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum | | For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung |
| | LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Eviter toute exposition au faisceau / Radiación Laser: Evite Exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden | | HOT SURFACE: Keep hands Away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden |
| | CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement Laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung | IPX-0 | CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen |
| | CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr | | |

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.