# CONSORT

---

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

**E802**

**E815**

**E831**

**E832**

**E833**

**E835**

**E844**

**E861**

**E862**

**E863**

**E865**

# QUALITYMASTERS

## CERTIFICATE OF REGISTRATION

This is to certify that:

### Consort NV

Parklaan 36  
B-2300 Turnhout

Has been assessed by QualityMasters in respect of its Quality Management System and found to comply with:

NEN-EN-ISO 9001 : 2000

Approval is hereby granted for registration providing the rules and conditions relating to certification are observed at all times.

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Lack of fulfillment of conditions as set forth in the certification regulations may render this certificate invalid. The use of the accreditation mark indicates accreditation in respect to the activities covered by the accreditation number EU9706010F.

# Q

# M

# REGISTRATION

# CERTIFICATE OF

Authorized signatory

*A. LEUWATY*  
QualityMasters B.V.

This certificate remains the property of QualityMasters B.V.

# WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory.

**CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse.

A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

## DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

**Electrophoresis Power Supply**

**content of the type numbers**

**E835, E831, E844, E802, E865, E861, E863, E815, E832, E833, E862**

to which this declaration relates is in conformity

with the following standards

**EN61010**

**LOW VOLTAGE DIRECTIVE 73/23/EEG**

**EN50081-1**

**EN50082-1**

**EN60555-2**

**EMC DIRECTIVE 89/336/EEG**

Turnhout, August 12, 1996

A handwritten signature in black ink, appearing to be 'J. G. G.', is written over a horizontal line.

on behalf of **CONSORT nv**

**Table of contents**

Introduction .....	5
Manufacturer .....	5
Warranty .....	5
Servicing.....	5
Specifications.....	6
Alarms .....	6
Very important .....	7
Pre-set limits .....	7
Keyboard .....	7
Programming.....	8

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**Forms**

Frequently Asked Questions.....	21
Technical assistance form .....	22
Your satisfaction .....	23

This instrument is manufactured with the latest technology and needs no particular maintenance. **CONSORT** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 4° and 40° C. Never store the instrument in a room with high humidity or at very low temperatures (condensation water!).

Connect the instrument only to an earthed power line. The required power source is indicated on the label at the back of the instrument. Do not cut and splice the power cord. When removing the power cord from the wall outlet, be sure to unplug by holding the plug attachment and not by pulling the cord. Do not hold the plug by wet hand.

---

**Manufacturer****CONSORT nv**

Parklaan 36  
B2300 Turnhout  
Belgium

Tel (+32)(14)41 12 79  
Fax (+32)(14)42 91 79  
E-mail: info@consort.be

---

**Warranty**

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **CONSORT** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **CONSORT** before returning any product for warranty repair on a freight prepaid basis!

**CONSORT** is not liable for consequential damages arising out of the use or handling of its products.

---

**Servicing**

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual  
Cables  
Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local **CONSORT** dealer.

<b>Ranges</b>	E835:	0...300 V, 0...500 mA, 0...150 W
	E831:	0...300 V, 0...1000 mA, 0...150 W
	E802:	0...300 V, 0...2000 mA, 0...300 W
	E844:	0...400 V, 0...400 mA, max. 80 W
	E863:	0...600 V, 0...250 mA, 0...150 W
	E865:	0...600 V, 0...500 mA, 0...150 W
	E861:	0...600 V, 0...1000 mA, 0...300 W
	E815:	0...1200 V, 0...500 mA, 0...300 W
	E832:	0...3000 V, 0...150 mA, 0...150 W
	E833:	0...3000 V, 0...300 mA, 0...300 W
E862:	0...6000 V, 0...150 mA, 0...300 W	

**Timer** 0:01...99:59 h  
0.01...99.99 kWh (not on E844)

**Programs** 9 programs with 9 steps

**Outputs** four, in parallel

**Leakage detection** > 500  $\mu$ A

**Minimum load resistance**

E835:	>20 $\Omega$
E831:	>10 $\Omega$
E844:	>30 $\Omega$
E802:	>5 $\Omega$
E863:	>70 $\Omega$
E865:	>30 $\Omega$
E861:	>15 $\Omega$
E815:	>70 $\Omega$
E832:	>600 $\Omega$
E833:	>300 $\Omega$
E862:	>1200 $\Omega$

**Ambient temp.** 4...40°C

**Relative humidity** 0...90 % (non-condensing!)

**Power supply** 210-250 V~ (105-125 V~), 50/60 Hz

**Cabinet** rugged plastic/metal cabinet

**Dimensions** 30 x 12 x 18 cm (80 and 150 W models)  
30 x 12 x 28 cm (300 W models and E832)

**Weight** 5 kg (80 and 150 W models)  
10 kg (300 W models)

### Alarms

**[ALrM][Gnd]** = A dangerous ground leakage of >500  $\mu$ A is present (check thoroughly the complete assembly).

**[ALrM][Shrt]** = Short-circuit condition! Connected tank has a too low resistance (check thoroughly the complete assembly).

This instrument is equipped with safety outputs on which the high voltage cannot suddenly appear. The built in microcomputer will always smoothly increase the voltage till one of the pre-set limits is reached. The unit is fully protected against any overload condition by a special safety system which automatically disconnects the AC line from the high voltage transformer. A flashing error code in the display and an acoustic alarm is activated whenever a dangerous condition is present.

This electrophoresis power supply is a high technology instrument available in several versions. As it is capable of giving dangerous voltage levels by which high power is involved, we suggest that you take a few moments to read this manual thoroughly. Although this instrument is equipped with all necessary safety features against abuse and other accidental failures, caution should be exercised when working with high voltage equipment. Therefore, **avoid to touch the outlets with any conducting object** and make sure there is a second person present for your safety in case of any severe electric shock. Never touch any part of the assembly (power supply, leads or tank) before having switched off. **Never manipulate with wet hands. Do not ground any of the outputs or the buffer in the tank.** Connect the outlets only to an insulated electrophoresis tank with safety cover. Never make any other connections, such as e.g. putting several power supplies in series or in parallel. In order to prevent electric shock, never open the back plate nor remove the cover. Do not expose the unit to rain or any other liquid. **Do not spill liquid or insert metal objects inside the unit.** Take care so that the power supply is not dropped to avoid damaging the cabinet which defeats safeguards or injuring yourself. If the unit has been dropped or the cabinet has been damaged, unplug it and have it checked by an authorised service technician to restore the safeguards. The fact that the unit operates satisfactorily does not imply that the unit is properly earthed or that it is completely safe. If in any doubt about the effective earthing of the unit, contact a qualified electrician. **Never block the ventilation holes** or place the unit in any enclosure unless proper ventilation is provided. Never place the unit near or over a radiator, heat register or stove. Avoid locations where the instrument is exposed directly to the sun light.

**Pre-set limits**

Parameters can be pre-set within the following limits:

E835:	3 ... 300 V,	5 ... 500 mA,	2 ... 150 W
E831:	3 ... 300 V,	10 ... 1000 mA,	2 ... 150 W
E802:	3 ... 300 V,	20 ... 2000 mA,	3 ... 300 W
E844:	4 ... 400 V,	4 ... 400 mA,	max. 80 W
E863:	6 ... 600 V,	3 ... 250 mA,	2 ... 150 W
E865:	6 ... 600 V,	5 ... 500 mA,	2 ... 150 W
E861:	6 ... 600 V,	10 ... 1000 mA,	3 ... 300 W
E815:	12 ... 1200 V,	5 ... 500 mA,	3 ... 300 W
E832:	30 ... 3000 V,	2 ... 150 mA,	2 ... 150 W
E833:	30 ... 3000 V,	3 ... 300 mA,	3 ... 300 W
E862:	60 ... 6000 V,	2 ... 150 mA,	3 ... 300 W

**Keyboard**

- ▲/▼** = Button for entering a value or for selecting a function.
- SET** = Button for programming the desired parameters, while the instrument is in stand-by.
- RUN** = Button to start an experiment.
- PROG** = Button to select a program.

1. Verify if the instrument is switched off.
  2. Connect the electrophoresis tank to the instrument.
  3. Switch the instrument on. If the timer is still running, the power supply will automatically proceed with the experiment for the remaining time. Press **SET** if you wish to program new parameters or want to go in stand-by during the run.
  4. The display shows [OFF]. Press **PROG** to select one of the 9 possible programs.
  5. The display shows e.g. [Px-y][ProG]. Select the desired program (x) and press **PROG**.
  6. The display shows e.g. [Px-y][ProG]. Select the desired step (y) and press **SET**.
  7. The display shows e.g. [3000][Px-y], while the constant voltage indicator is blinking. Enter the desired voltage (V) limit and press **SET**.
  8. The display shows e.g. [Px-y][150], while the constant current indicator is blinking. Enter the desired current (mA) limit and press **SET**.
  9. The display shows e.g. [Px-y][300], while the constant power indicator is blinking (not on E844). Enter the desired power (W) limit and press **SET**.
  10. The display shows e.g. [Px-y][hour] (not on E844). Select the desired time-units and press **SET**.
  11. The display shows e.g. [Px-y][99.59], while the time-indicator is blinking. Enter the desired time and press **SET**. Pre-set to [--.--] if no timing is required.
  12. The instrument will let you choose between [Px-y][End] or [Px-y][Cnt]. The experiment will automatically proceed with a next set of parameters (multiple step programming), when [Cnt] is selected. [End] means that the experiment will stop after this step. Select the desired procedure and press **SET**.
  13. Press **RUN** to start the experiment.
- Press  $\blacktriangle$  or  $\blacktriangledown$  to display (right) the current, power, time, or program number.
  - At any time, you can avoid programming unnecessary parameters. To do so, simply press **RUN** instead of **SET** after entering the last relevant value and the instrument will immediately start the experiment keeping the other parameters to their previous values.
  - Important formulae determine which parameter will be kept constant:  
 Power (W) = Voltage (V) x Current (A)  
 Resistance ( $\Omega$ ) = Voltage (V) / Current (A)
  - To ensure that the desired parameter is kept constant, we suggest to pre-set all other parameters to maximum.
  - The automatic crossover indicators will only illuminate (blinking units) when the regulation is complete and one of the parameters is kept constant.
  - Press **PROG** during a run to see the pre-set parameters of the displayed values.
  - Press **RUN** during a run to see the total elapsed time of the experiment.

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**What are the relations between Voltage, Current, Power and Resistance?**

Power = Voltage x Current (1)

Resistance = Voltage / Current (2)

**How does a power supply react after pressing RUN?**

The internal generator will start building up the high voltage at the output terminals while voltage and current are constantly measured and power calculated. When one of the pre-set parameters is exceeded, the generator stops and will keep that parameter constant.

**How important is the resistance of an electrophoresis unit?**

The resistance of an electrophoresis unit depends its size, gel thickness, amount of buffer, buffer conductivity and temperature. This resistance will normally decrease in time due to a slowly increasing temperature. Electrophoresis units which have a resistance below the minimum load resistance of a power supply will trigger an alarm! Read the output voltage and current during a run to measure the resistance and use formula (2) to calculate the value.

**How to keep a constant voltage during a run?**

Program the desired voltage and a higher current and power than the maximum expected values:

Current > Voltage / Resistance

Power > Voltage x Current

**How to keep a constant current during a run?**

Program the desired current and a higher voltage and power than the maximum expected values:

Voltage > Current x Resistance

Power > Voltage x Current

**How to keep a constant power during a run?**

Program the desired power and a higher voltage and current than the maximum expected values:

Voltage > Current x Resistance

Current > Voltage / Resistance

**Why are my output values different from those of a similar experiment?**

Either your programmed parameters are not equal to those described or the resistance of your electrophoresis unit is different (see above). It cannot be due to e.g. an other model of power supply as the relations between Voltage, Current, Power and Resistance are monitored in the same way by any instrument (the electrical laws cannot be disregarded!).

**What about connecting more than one unit to the same power supply?**

The outlets being in parallel each electrophoresis unit will be supplied with exactly the same voltage. However, current and power may differ due to differences between them even when exactly the same model, gel, buffers, etc... are used. Therefore, it is recommended to run several electrophoresis units on the same power supply only in the constant voltage mode.

If you need assistance with a technical problem, make sure to fill the **complete** checklist (26 points!) and fax a copy to Consort: +32/14/429179.

1. purchase date:  2. purchased from:

2. indicate model:  4. serial no:

5. organisation:

6. name:

7. division:

8. address:

9. postal code:  10. city:

11. telephone:  12. telefax:

13. e-mail:

14. type of connected electrophoresis unit:

**Perform the following steps:**

15. read the preset voltage.	V:	<input type="text"/>
16. read the preset current.	mA:	<input type="text"/>
17. read the preset power.	W:	<input type="text"/>
18. connect the electrophoresis unit (or a dummy resistance) and press RUN.	yes?	<input type="text"/>
19. wait until one of the units blinks, then read the output voltage.	V:	<input type="text"/>
20. read the output current.	mA:	<input type="text"/>
21. read the output power.	W:	<input type="text"/>
22. reset the power supply by switching on while holding SET pressed.	yes?	<input type="text"/>
23. leave ALL outputs open (no load connected).	yes?	<input type="text"/>
24. program the voltage to maximum.	V:	<input type="text"/>
25. press RUN and wait until the V-unit blinks, then read the output voltage.	V:	<input type="text"/>
26. if any alarm occurs, what is shown on both displays?	<input type="text"/>	<input type="text"/>