




## MAGNETIC SEPARATION

BIOVENDIS

### M 12+12



Consists of a clear polyacril block with 12× extra strong permanent magnets for a fast magnetic separation with 1.5 ml and 2.0 ml tubes.

- Optimized for Kits
- With 48 tube positions including 24 positions for magnetic separation (12× for 1.5 ml and 12× for 2.0 ml tubes)
- For 1.5 ml (conical) and 2.0 ml (conical recommended) tubes
- Transparent material allows constant visual control of the separation
- For fast separation (< 2 min)

**IMPORTANT!** We recommend cleaning the Magnetic Stand M96 with deionized water only; alcohol-based cleaning agents may cause cracking of the stand material.

**For Research Use Only. Not for human or animal therapeutic or diagnostic use.**

	<b>M 12+12</b>
Cat. No.	<b>210141</b>
Format	1.5 / 2.0ml tubes
Plastic	clear polyacril
Colour	<input type="checkbox"/>
Items per Unit	1
Criteria	for 12× 1.5 and 2.0ml tubes